







Terufusion®







Infusion System

Infusion and Extension Sets






Infusion Sets



Type	General					
Product code	TI*PU300WY	TI*PJ300WY01	YI*PU400WN01	YI*PU400WY51	YI*PU400WY21	TI*PM370WYE ¹
AFF clip	Yes	Yes	No	Yes	Yes	Yes
Filter (um)	15	15	15	15	15	15
Pump tube (mm)	2000	2000	2000	2500	2200	2200
Injector port	None	Needle port	Needle port	Needle port	Needle port	NLCD ²
Drip rate (drops/ml)	20	20	20	20	20	60
Material - tubing	PVC ³	Polybutadiene ⁴	PVC ³	PVC ³	PVC ³	PVC ³
Packaging (unit/carton)	20/200	20/200	20/200	20/200	20/200	20/200
Priming volume (ml)	15,7	16,8	16,9	20,7	18,1	19,4
Colour ⁵	Transparent	Transparent	Transparent	Transparent	Transparent	Transparent
Pressure resistance (bar)	2	2	2	2	2	2
Total length of sets	2150	2440	2270	2770	2470	2985
						

Type	General	General (amber)		Oncology	Transfusion	
Product code	SP*PT306WY01E	TI*PA300WN01E	TI*PA300WY01E	SP*PF306WY01E	TB*PU300SE	TB*PU300SYE
AFF clip	Yes	No	Yes	Yes	No	Yes
Filter (um)	15	15	15	0,2	200	200
Tube length (mm)	2200	2000	2000	2200	2200	2200
Injector port	NLCD ²	Needle port	Needle port	NLCD ²	None	None
Drip rate (drops/ml)	20	20	20	20	20	20
Material - tubing	Polybutadiene ⁴	PVC ³	PVC ³	Polybutadiene ⁴	PVC ³	PVC ³
Packaging (unit/carton)	20/200	20/200	20/200	20/200	20/200	20/200
Priming volume (ml)	18,3	17,4	17,4	20,2	17,4	17,4
Colour ⁵	Transparent	Amber	Amber	Transparent	Transparent	Transparent
Pressure resistance (bar)	2	2	2	2	2	2
Total length of sets	2640	2640	2640	2860	2420	2420
						

Extension Sets

Type	Basic					
Product code	EX+300TN1FL	EX+200TN1FL	EX+150TN1FL	EX+075TN1FL	EX+050TN1FL	EX+025TN1FL
Material - tubing	PVC ¹	PVC ¹	PVC ¹	PVC ¹	PVC ¹	PVC ¹
Colour - tubing	Transparent	Transparent	Transparent	Transparent	Transparent	Transparent
Tube length (cm)	300	200	150	75	50	25
Connection 1	FFLL ²	FFLL ²	FFLL ²	FFLL ²	FFLL ²	FFLL ²
Connection 2	FMLL ³	FMLL ³	FMLL ³	FMLL ³	FMLL ³	FMLL ³
Clip	No	No	No	No	No	No
Priming volume (ml)	5,3	3,5	2,6	1,5	1	0,5
Pressure resistance (bar)	4	4	4	4	4	4
						

Type	Amber				3-way tap			
Product code	EX+250AN1FL	EX+200AN1FL	EX+150AN1FL	EX+025AN1FL	EX+200TNTAP	EX+150TNTAP	EX+050TNTAP	EX+025TNTAP
Material - tubing	PUR ⁴	PUR ⁴	PUR ⁴	PUR ⁴	PVC ¹	PVC ¹	PVC ¹	PVC ¹
Colour - tubing	Amber	Amber	Amber	Amber	Transparent	Transparent	Transparent	Transparent
Tube length (cm)	250	200	150	25	200	150	50	25
Connection 1	FFLL ²	FFLL ²	FFLL ²	FFLL ²	3-way tap	3-way tap	3-way tap	3-way tap
Connection 2	RMLL ⁵	RMLL ⁵	RMLL ⁵	RMLL ⁵	FMLL ³	FMLL ³	FMLL ³	FMLL ³
Clip	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Priming volume (ml)	2,0	1,6	1,2	0,5	2,2	2,0	1,0	0,5
Pressure resistance (bar)	21	21	21	21	4	4	4	4
								

Type	Add-on	
Product code	A0+T2NFD	A0+T4NFD
Material - tubing	PVC ¹	PVC ¹
Colour - tubing	Transparent	Transparent
Tube length (cm)	25+3	25+3
Connection 1	2 x (NLCD ⁶ +BCV ⁷)	4 x (NLCD ⁶ +BCV ⁷)
Connection 2	Piercing connector	Piercing connector
Clip	Yes	Yes
Priming volume (ml)	2,0	2,4
Pressure resistance (bar)	2	2
		



Terumo Infusion Sets

Easy handling

- Easy setting following the instruction screen of the pump
- Sharp 2-way spike penetrating bags and bottles easily
- Adjustable roller clamp with notch to mount tubing during priming


Increased patient safety

- AFF clip for prevention of accidental bolus infusion
- Does not contain DEHP
- Y-site with NLCD* for prevention of needle stick injuries

Highly accurate infusion

- Reliable high accuracy combination of Terumo pump & sets**
- Durable tubes for continuous accurate infusion

* Needleless connecting device
 ** Flow rate accuracy: within ± 5%



Full focus infusion system

Terumo offers a comprehensive portfolio of solutions for the complete infusion line-up: user-friendly syringe and volumetric pumps and racks, as well as a wide range of highly accurate, easy-to-use disposables, such as an extended range of infusion sets and syringes.



Terumo Corporation
 +81 3 3374 8111

Terumo Europe NV
 +32 16 38 12 11

Your Terumo partner